

Restores Skin's Elasticity, Reduces Wrinkles, and Tones the Facial Muscles

With GEM Pro you can restore skin's vitality and flexibility, reduce fine lines, tone and improve skin texture. Portable, lightweight, and easy to operate - this beauty product is a revolutionary new system for Spa, Salon, and Home use.

GEM Pro comprises 3 beneficial skin care programs using 3 types of electric currents:

- 1. GALVANIC CURRENT:** A continuous current which increases the absorption of skin care products and improves skin appearance.
- 2. EMS (ELECTRICAL MUSCLE STIMULATION):** Improves vitality and tones facial muscles.
- 3. MICRO-CURRENT:** Variable alternate current that minimizes the appearance of wrinkles.



Function 1: Seven Golden Needles Microcurrent

Thermal microcurrent, as an anti-wrinkle treatment, is applied with 7 small round-tip contact points, that place electrodes on the deepest zones where wrinkles are formed. This thermal microcurrent stimulation produces proliferation of fibroblasts cells on the skin. Fibroblasts help to generate collagen and elastin. After the anti-wrinkle treatment deep wrinkles will form mild swelling of skin elastic fibers that will promote toning. The depth of wrinkles will rise and restore its original smooth appearance. With further treatment, efficiency becomes apparent and wrinkles form a new line of soft collagen. This treatment can also be used to prevent the appearance of wrinkles.



Function 2: Germanium Probe Electrode

The facial muscle toning treatment consists of applying microcurrent on flaccid skin of the lower face and neck (neck, cheeks, and top of the chin). The treatment acts directly on skin cell function, muscles, and lymphatic system.

Use of electric stimulation with low-frequency micro current is used to:

- Restore each layer of skin
- Regulate skin microcirculation and lymphatic drainage
- Stimulate collagen and elastin production
- Stimulate facial muscles to eliminate excess fat deposits on face and neck, double chin.



Function 2: Germanium Y-Shape Applicator

LIFTING THERAPY (Lifting Ion +/-)

Continuous current (CC) application helps the material exchange process within the skin to improve its vitality, control oil secretion, enrich skin nutrition, moisturize and regenerate, and to stimulate special production of collagen and elastin within fibroblasts. Collagen forms the internal framework of the skin. It strengthens its elasticity, vitality and moistening capabilities. The elastic fibers improve flexibility and restore skin plasticity. Lifting therapy or iontophoresis consists of penetrating skin care products into the skin cells. There are two types of lifting therapy:



Negative Ion Lifting / DEEP CLEANSING
(Y shape applicator with - electrode):

Positive Ion Lifting / DEEP MOISTURIZING
(Y shape applicator with + electrode):

